

Antibodies Product Datasheet

Product Name:	DNA pol ϵ B Antibody
Catalog Number:	C15410
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human DNA pol ϵ B.
Concentration:	1 mg/ml
Gene Name:	POLE2
Protein Name:	DNA polymerase epsilon subunit 2
Other Name:	POLE2; DPE2; DNA polymerase epsilon subunit 2; DNA polymerase II subunit 2; DNA polymerase epsilon subunit B
Human Swiss Prot No.:	P56282
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DNA pol ϵ B Polyclonal Antibody detects endogenous levels of DNA pol ϵ B protein.
Storage and Stability:	-20°C/1 year