

Antibodies Product Datasheet

Product Name:	DNA pol δ 3 Antibody
Catalog Number:	C15428
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human POLD3. AA range:411-460
Concentration:	1 mg/ml
Gene Name:	POLD3
Protein Name:	DNA polymerase delta subunit 3
Other Name:	POLD3; KIAA0039; DNA polymerase delta subunit 3; DNA polymerase delta subunit p66
Human Swiss Prot No.:	Q15054
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DNA pol δ 3 Polyclonal Antibody detects endogenous levels of DNA pol δ 3 protein.
Storage and Stability:	-20°C/1 year