

Antibodies Product Datasheet

Product Name:	CD3EAP Antibody
Catalog Number:	C15463
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD3EAP. AA range:441-490
Concentration:	1 mg/ml
Gene Name:	CD3EAP
Protein Name:	DNA-directed RNA polymerase I subunit RPA34
Other Name:	CD3EAP; ASE1; CAST; PAF49; DNA-directed RNA polymerase I subunit RPA34; A34.5; Antisense to ERCC-1 protein; ASE-1; CD3-epsilon-associated protein; CAST; CD3E-associated protein; RNA polymerase I-associated factor PAF49
Human Swiss Prot No.:	O15446
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CD3EAP Polyclonal Antibody detects endogenous levels of CD3EAP protein.
Storage and Stability:	-20°C/1 year