

Antibodies Product Datasheet

Product Name:	POLR2J1 Antibody
Catalog Number:	C15467
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPB11. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	POLR2J
Protein Name:	DNA-directed RNA polymerase II subunit RPB11-a
Other Name:	POLR2J; POLR2J1; DNA-directed RNA polymerase II subunit RPB11-a; RNA polymerase II subunit B11-a; RPB11a; DNA-directed RNA polymerase II subunit J-1; RNA polymerase II 13.3 kDa subunit P52435
Human Swiss Prot No.:	P52435
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	POLR2J1 Polyclonal Antibody detects endogenous levels of POLR2J1 protein.
Storage and Stability:	-20°C/1 year