

Antibodies Product Datasheet

Product Name:	POLR3D Antibody
Catalog Number:	C15477
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RPC4. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	POLR3D
Protein Name:	DNA-directed RNA polymerase III subunit RPC4
Other Name:	POLR3D; BN51; BN51T; DNA-directed RNA polymerase III subunit RPC4; RNA polymerase III subunit C4; DNA-directed RNA polymerase III subunit D; Protein BN51; RNA polymerase III 47 kDa subunit; RPC53 homolog
Human Swiss Prot No.:	P05423
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	POLR3D Polyclonal Antibody detects endogenous levels of POLR3D protein.
Storage and Stability:	-20°C/1 year