

Antibodies Product Datasheet

Product Name:	CSP Antibody
Catalog Number:	C15493
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DNAJC5. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	DNAJC5
Protein Name:	DnaJ homolog subfamily C member 5
Other Name:	DNAJC5; CSP; DnaJ homolog subfamily C member 5; Cysteine string protein; CSP
Human Swiss Prot No.:	Q9H3Z4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CSP Polyclonal Antibody detects endogenous levels of CSP protein.
Storage and Stability:	-20°C/1 year