

Antibodies Product Datasheet

Product Name:	Dynein IC1 Antibody
Catalog Number:	C15510
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DNAI1. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	DNAI1
Protein Name:	Dynein intermediate chain 1 axonemal
Other Name:	DNAI1; Dynein intermediate chain 1; axonemal; Axonemal dynein intermediate chain 1
Human Swiss Prot No.:	Q9UI46
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dynein IC1 Polyclonal Antibody detects endogenous levels of Dynein IC1 protein.
Storage and Stability:	-20°C/1 year