

Antibodies Product Datasheet

Product Name:	Dynein LC 2B Antibody
Catalog Number:	C15517
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DYNLRB2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	DYNLRB2
Protein Name:	Dynein light chain roadblock-type 2
Other Name:	DYNLRB2; DNCL2B; DNLC2B; ROBLD2; Dynein light chain roadblock-type 2; Dynein light chain 2B; cytoplasmic; Roadblock domain-containing protein 2
Human Swiss Prot No.:	Q8TF09
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dynein LC 2B Polyclonal Antibody detects endogenous levels of Dynein LC 2B protein.
Storage and Stability:	-20°C/1 year