

Antibodies Product Datasheet

Product Name:	HEI10 Antibody
Catalog Number:	C15533
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCNB1IP1. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	CCNB1IP1
Protein Name:	E3 ubiquitin-protein ligase CCNB1IP1
Other Name:	CCNB1IP1; C14orf18; HEI10; E3 ubiquitin-protein ligase CCNB1IP1; Cyclin-B1-interacting protein 1; Human enhancer of invasion 10
Human Swiss Prot No.:	Q9NPC3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HEI10 Polyclonal Antibody detects endogenous levels of HEI10 protein.
Storage and Stability:	-20°C/1 year