

Antibodies Product Datasheet

Product Name:	MARCH1 Antibody
Catalog Number:	C15546
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human 40603.
Concentration:	1 mg/ml
Gene Name:	MARCH1
Protein Name:	E3 ubiquitin-protein ligase MARCH1
Other Name:	MARCH1; RNF171; E3 ubiquitin-protein ligase MARCH1; Membrane-associated RING finger protein 1; Membrane-associated RING-CH protein 1; MARCH-I; RING finger protein 171
Human Swiss Prot No.:	Q8TCQ1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MARCH1 Polyclonal Antibody detects endogenous levels of 40603 protein.
Storage and Stability:	-20°C/1 year