

Antibodies Product Datasheet

Product Name:	MARCH4 Antibody
Catalog Number:	C15549
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MARCH4. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	MARCH4
Protein Name:	E3 ubiquitin-protein ligase MARCH4
Other Name:	MARCH4; KIAA1399; RNF174; E3 ubiquitin-protein ligase MARCH4; Membrane-associated RING finger protein 4; Membrane-associated RING-CH protein IV; MARCH-IV; RING finger protein 174
Human Swiss Prot No.:	Q9P2E8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MARCH4 Polyclonal Antibody detects endogenous levels of 40606 protein.
Storage and Stability:	-20°C/1 year