

Antibodies Product Datasheet

Product Name:	MYLIP Antibody
Catalog Number:	C15558
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYLIP. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	MYLIP
Protein Name:	E3 ubiquitin-protein ligase MYLIP
Other Name:	MYLIP; BZF1; IDOL; BM-023; PP5242; E3 ubiquitin-protein ligase MYLIP; Inducible degrader of the LDL-receptor; Idol; Myosin regulatory light chain interacting protein; MIR
Human Swiss Prot No.:	Q8WY64
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MYLIP Polyclonal Antibody detects endogenous levels of MYLIP protein.
Storage and Stability:	-20°C/1 year