

Antibodies Product Datasheet

Product Name:	RNF138 Antibody
Catalog Number:	C15572
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RNF138. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	RNF138
Protein Name:	E3 ubiquitin-protein ligase RNF138
Other Name:	RNF138; NARF; HSD-4; HSD4; E3 ubiquitin-protein ligase RNF138; Nemo-like kinase-associated RING finger protein; NLK-associated RING finger protein; hNARF; RING finger protein 138
Human Swiss Prot No.:	Q8WVD3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RNF138 Polyclonal Antibody detects endogenous levels of RNF138 protein.
Storage and Stability:	-20°C/1 year