

Antibodies Product Datasheet

Product Name:	EF-1 γ Antibody
Catalog Number:	C15610
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EEF1G. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	EEF1G
Protein Name:	Elongation factor 1-gamma
Other Name:	EEF1G; EF1G; Elongation factor 1-gamma; EF-1-gamma; eEF-1B gamma
Human Swiss Prot No.:	P26641
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EF-1 γ Polyclonal Antibody detects endogenous levels of EF-1 γ protein.
Storage and Stability:	-20°C/1 year