

Antibodies Product Datasheet

Product Name:	ELOVL6 Antibody
Catalog Number:	C15622
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELOVL6. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	ELOVL6
Protein Name:	Elongation of very long chain fatty acids protein 6
Other Name:	ELOVL6; FACE; LCE; Elongation of very long chain fatty acids protein 6; 3-keto acyl-CoA synthase ELOVL6; ELOVL fatty acid elongase 6; ELOVL FA elongase 6; Fatty acid elongase 2; hELO2; Fatty acyl-CoA elongase; Long-chain fatty-acyl elongase
Human Swiss Prot No.:	Q9H5J4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ELOVL6 Polyclonal Antibody detects endogenous levels of ELOVL6 protein.
Storage and Stability:	-20°C/1 year