

Antibodies Product Datasheet

Product Name:	Endo G-L1 Antibody
Catalog Number:	C15631
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ENDOGL1. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	EXOG
Protein Name:	Nuclease EXOG mitochondrial
Other Name:	EXOG; ENDOGL1; ENDOGL2; ENGL; Nuclease EXOG; mitochondrial; Endonuclease G-like 1; Endo G-like 1
Human Swiss Prot No.:	Q9Y2C4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Endo G-L1 Polyclonal Antibody detects endogenous levels of Endo G-L1 protein.
Storage and Stability:	-20°C/1 year