

Antibodies Product Datasheet

Product Name:	MAN1B1 Antibody
Catalog Number:	C15636
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MAN1B1 . at AA range: 100-180
Concentration:	1 mg/ml
Gene Name:	MAN1B1
Protein Name:	Endoplasmic reticulum mannosyl-oligosaccharide 1,2-alpha-mannosidase
Other Name:	MAN1B1; Endoplasmic reticulum mannosyl-oligosaccharide 1; 2-alpha-mannosidase; ER alpha-1,2-mannosidase; ER mannosidase 1; ERMan1; Man9GlcNAc2-specific-processing alpha-mannosidase; Mannosidase alpha class 1B member 1
Human Swiss Prot No.:	Q9UKM7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAN1B1 Polyclonal Antibody detects endogenous levels of MAN1B1 protein.
Storage and Stability:	-20°C/1 year