

Antibodies Product Datasheet

Product Name:	Epsin 1 Antibody
Catalog Number:	C15666
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Epsin 1 . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	EPN1
Protein Name:	Epsin-1
Other Name:	EPN1; Epsin-1; EH domain-binding mitotic phosphoprotein; EPS-15-interacting protein 1
Human Swiss Prot No.:	Q9Y6I3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Epsin 1 Polyclonal Antibody detects endogenous levels of Epsin 1 protein.
Storage and Stability:	-20°C/1 year