

Antibodies Product Datasheet

Product Name:	eIF1B Antibody
Catalog Number:	C15693
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human eIF1B.
Concentration:	1 mg/ml
Gene Name:	EIF1B
Protein Name:	Eukaryotic translation initiation factor 1b
Other Name:	EIF1B; Eukaryotic translation initiation factor 1b; eIF1b; Protein translation factor SUI1 homolog GC20
Human Swiss Prot No.:	O60739
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF1B Polyclonal Antibody detects endogenous levels of eIF1B protein.
Storage and Stability:	-20°C/1 year