

Antibodies Product Datasheet

Product Name:	eIF3 α Antibody
Catalog Number:	C15705
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from eIF3 α . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	EIF3J
Protein Name:	Eukaryotic translation initiation factor 3 subunit J
Other Name:	EIF3J; EIF3S1; Eukaryotic translation initiation factor 3 subunit J; eIF3j; Eukaryotic translation initiation factor 3 subunit 1; eIF-3-alpha; eIF3 p35
Human Swiss Prot No.:	O75822
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF3 α Polyclonal Antibody detects endogenous levels of eIF3 α protein.
Storage and Stability:	-20°C/1 year