

Antibodies Product Datasheet

Product Name:	eIF3ε Antibody
Catalog Number:	C15711
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EIF3F. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	EIF3F
Protein Name:	Eukaryotic translation initiation factor 3 subunit F
Other Name:	EIF3F; EIF3S5; Eukaryotic translation initiation factor 3 subunit F; eIF3f; Deubiquitinating enzyme eIF3f; Eukaryotic translation initiation factor 3 subunit 5; eIF-3-epsilon; eIF3 p47
Human Swiss Prot No.:	O00303
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF3ε Polyclonal Antibody detects endogenous levels of eIF3ε protein.
Storage and Stability:	-20°C/1 year