

Antibodies Product Datasheet

Product Name:	eIF3 ζ Antibody
Catalog Number:	C15714
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EIF3D. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	EIF3D
Protein Name:	Eukaryotic translation initiation factor 3 subunit D
Other Name:	EIF3D; EIF3S7; Eukaryotic translation initiation factor 3 subunit D; eIF3d; Eukaryotic translation initiation factor 3 subunit 7; eIF-3-zeta; eIF3 p66
Human Swiss Prot No.:	O15371
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF3 ζ Polyclonal Antibody detects endogenous levels of eIF3 ζ protein.
Storage and Stability:	-20°C/1 year