

Antibodies Product Datasheet

Product Name:	eIF4E3 Antibody
Catalog Number:	C15723
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EIF4E3. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	EIF4E3
Protein Name:	Eukaryotic translation initiation factor 4E type 3
Other Name:	EIF4E3; Eukaryotic translation initiation factor 4E type 3; eIF-4E type 3; eIF-4E3; eIF4E type 3; eIF4E-3
Human Swiss Prot No.:	Q8N5X7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF4E3 Polyclonal Antibody detects endogenous levels of eIF4E3 protein.
Storage and Stability:	-20°C/1 year