

Antibodies Product Datasheet

Product Name:	eIF5 Antibody
Catalog Number:	C15725
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human eIF5.
Concentration:	1 mg/ml
Gene Name:	EIF5
Protein Name:	Eukaryotic translation initiation factor 5
Other Name:	EIF5; Eukaryotic translation initiation factor 5; eIF-5
Human Swiss Prot No.:	P55010
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF5 Polyclonal Antibody detects endogenous levels of eIF5 protein.
Storage and Stability:	-20°C/1 year