

Antibodies Product Datasheet

Product Name:	eIF5B Antibody
Catalog Number:	C15728
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EIF5B. AA range:1051-1100
Concentration:	1 mg/ml
Gene Name:	EIF5B
Protein Name:	Eukaryotic translation initiation factor 5B
Other Name:	EIF5B; IF2; KIAA0741; Eukaryotic translation initiation factor 5B; eIF-5B; Translation initiation factor IF-2
Human Swiss Prot No.:	O60841
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF5B Polyclonal Antibody detects endogenous levels of eIF5B protein.
Storage and Stability:	-20°C/1 year