

Antibodies Product Datasheet

Product Name:	FBP2 Antibody
Catalog Number:	C15742
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human FBP2.
Concentration:	1 mg/ml
Gene Name:	KHSRP
Protein Name:	Far upstream element-binding protein 2
Other Name:	KHSRP; FUBP2; Far upstream element-binding protein 2; FUSE-binding protein 2; KH type-splicing regulatory protein; KSRP; p75
Human Swiss Prot No.:	Q92945
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FBP2 Polyclonal Antibody detects endogenous levels of FBP2 protein.
Storage and Stability:	-20°C/1 year