

Antibodies Product Datasheet

Product Name:	FAIM2 Antibody
Catalog Number:	C15744
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human FAIM2.
Concentration:	1 mg/ml
Gene Name:	FAIM2
Protein Name:	Protein lifeguard 2
Other Name:	FAIM2; KIAA0950; LFG; LFG2; NMP35; TMBIM2; Protein lifeguard 2; Fas apoptotic inhibitory molecule 2; Neural membrane protein 35; Transmembrane BAX inhibitor motif-containing protein 2
Human Swiss Prot No.:	Q9BWQ8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAIM2 Polyclonal Antibody detects endogenous levels of FAIM2 protein.
Storage and Stability:	-20°C/1 year