

Antibodies Product Datasheet

Product Name:	GRIN1 Antibody
Catalog Number:	C15835
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPRIN1. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	GPRIN1
Protein Name:	G protein-regulated inducer of neurite outgrowth 1
Other Name:	GPRIN1; KIAA1893; G protein-regulated inducer of neurite outgrowth 1; GRIN1
Human Swiss Prot No.:	Q7Z2K8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRIN1 Polyclonal Antibody detects endogenous levels of GRIN1 protein.
Storage and Stability:	-20°C/1 year