

Antibodies Product Datasheet

Product Name:	GRIN3 Antibody
Catalog Number:	C15837
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPRIN3. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	GPRIN3
Protein Name:	G protein-regulated inducer of neurite outgrowth 3
Other Name:	GPRIN3; KIAA2027; G protein-regulated inducer of neurite outgrowth 3; GRIN3
Human Swiss Prot No.:	Q6ZVF9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRIN3 Polyclonal Antibody detects endogenous levels of GRIN3 protein.
Storage and Stability:	-20°C/1 year