

Antibodies Product Datasheet

Product Name:	FUT2 Antibody
Catalog Number:	C15839
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human FUT2.
Concentration:	1 mg/ml
Gene Name:	FUT2
Protein Name:	Galactoside 2-alpha-L-fucosyltransferase 2
Other Name:	FUT2; SEC2; Galactoside 2-alpha-L-fucosyltransferase 2; Alpha(1; 2)FT 2; Fucosyltransferase 2; GDP-L-fucose:beta-D-galactoside 2-alpha-L-fucosyltransferase 2; SE2; Secretor blood group alpha-2-fucosyltransferase; Secretor factor; Se
Human Swiss Prot No.:	Q10981
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FUT2 Polyclonal Antibody detects endogenous levels of FUT2 protein.
Storage and Stability:	-20°C/1 year