

Antibodies Product Datasheet

Product Name:	FucT-III Antibody
Catalog Number:	C15840
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FUT3. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	FUT3
Protein Name:	Galactoside 3(4)-L-fucosyltransferase
Other Name:	FUT3; FT3B; LE; Galactoside 3(4)-L-fucosyltransferase; Blood group Lewis alpha-4-fucosyltransferase; Lewis FT; Fucosyltransferase 3; Fucosyltransferase III; FucT-III
Human Swiss Prot No.:	P21217
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FucT-III Polyclonal Antibody detects endogenous levels of FucT-III protein.
Storage and Stability:	-20°C/1 year