

Antibodies Product Datasheet

Product Name:	GABAA R α 4 Antibody
Catalog Number:	C15850
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GABRA4. AA range:81-130
Concentration:	1 mg/ml
Gene Name:	GABRA4
Protein Name:	Gamma-aminobutyric acid receptor subunit alpha-4
Other Name:	GABRA4; Gamma-aminobutyric acid receptor subunit alpha-4; GABA(A) receptor subunit alpha-4
Human Swiss Prot No.:	P48169
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABAA R α 4 Polyclonal Antibody detects endogenous levels of GABAA R α 4 protein.
Storage and Stability:	-20°C/1 year