

Antibodies Product Datasheet

Product Name:	GABAA Rε Antibody
Catalog Number:	C15857
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human GABAA range: Rε.
Concentration:	1 mg/ml
Gene Name:	GABRE
Protein Name:	Gamma-aminobutyric acid receptor subunit epsilon
Other Name:	GABRE; Gamma-aminobutyric acid receptor subunit epsilon; GABA(A) receptor subunit epsilon
Human Swiss Prot No.:	P78334
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABAA Rε Polyclonal Antibody detects endogenous levels of GABAA Rε protein.
Storage and Stability:	-20°C/1 year