

## Antibodies Product Datasheet

|                               |                                                                                                                                       |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | IFI-16 Antibody                                                                                                                       |
| <b>Catalog Number:</b>        | C15867                                                                                                                                |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                     |
| <b>Applications:</b>          | WB, ELISA                                                                                                                             |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.                                                   |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                 |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human IFI16. AA range:731-780                                     |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                               |
| <b>Gene Name:</b>             | IFI16                                                                                                                                 |
| <b>Protein Name:</b>          | Gamma-interferon-inducible protein 16                                                                                                 |
| <b>Other Name:</b>            | IFI16; IFNGIP1; Gamma-interferon-inducible protein 16; Ifi-16; Interferon-inducible myeloid differentiation transcriptional activator |
| <b>Human Swiss Prot No.:</b>  | Q16666                                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                               |
| <b>Specificity:</b>           | IFI-16 Polyclonal Antibody detects endogenous levels of IFI-16 protein.                                                               |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                          |