

Antibodies Product Datasheet

Product Name:	GCP4 Antibody
Catalog Number:	C15872
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TUBGCP4. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	TUBGCP4
Protein Name:	Gamma-tubulin complex component 4
Other Name:	TUBGCP4; 76P; GCP4; Gamma-tubulin complex component 4; GCP-4; hGCP4; Gamma-ring complex protein 76 kDa; h76p; hGrip76
Human Swiss Prot No.:	Q9UGJ1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GCP4 Polyclonal Antibody detects endogenous levels of GCP4 protein.
Storage and Stability:	-20°C/1 year