

Antibodies Product Datasheet

Product Name:	Connexin 25 Antibody
Catalog Number:	C15889
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GJB7. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	GJB7
Protein Name:	Gap junction beta-7 protein
Other Name:	GJB7; CX25; Gap junction beta-7 protein; Connexin-25; Cx25
Human Swiss Prot No.:	Q6PEY0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Connexin 25 Polyclonal Antibody detects endogenous levels of Connexin 25 protein.
Storage and Stability:	-20°C/1 year