

Antibodies Product Datasheet

Product Name:	Connexin 31.3 Antibody
Catalog Number:	C15890
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GJC3. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	GJC3
Protein Name:	Gap junction gamma-3 protein
Other Name:	GJC3; GJE1; Gap junction gamma-3 protein; Connexin-30.2; Cx30.2; Connexin-31.3; Cx31.3; Gap junction epsilon-1 protein
Human Swiss Prot No.:	Q8NFK1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Connexin 31.3 Polyclonal Antibody detects endogenous levels of Connexin 31.3 protein.
Storage and Stability:	-20°C/1 year