

Antibodies Product Datasheet

Product Name:	GFR α -4 Antibody
Catalog Number:	C15896
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GFRA4. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	GFRA4
Protein Name:	GNDF family receptor alpha-4
Other Name:	GFRA4; GDNF family receptor alpha-4; GDNF receptor alpha-4; GDNFR-alpha-4; GFR-alpha-4; Persephin receptor
Human Swiss Prot No.:	Q9GZZ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GFR α -4 Polyclonal Antibody detects endogenous levels of GFR α -4 protein.
Storage and Stability:	-20°C/1 year