

Antibodies Product Datasheet

Product Name:	FX Antibody
Catalog Number:	C15900
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TSTA3. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	TSTA3
Protein Name:	GDP-L-fucose synthase
Other Name:	TSTA3; SDR4E1; GDP-L-fucose synthase; GDP-4-keto-6-deoxy-D-mannose-3; 5-epimerase-4-reductase; Protein FX; Red cell NADP(H)-binding protein; Short-chain dehydrogenase/reductase family 4E member 1
Human Swiss Prot No.:	Q13630
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FX Polyclonal Antibody detects endogenous levels of FX protein.
Storage and Stability:	-20°C/1 year