

Antibodies Product Datasheet

Product Name:	GSTT1/4 Antibody
Catalog Number:	C15930
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GSTT1/4. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	GSTT1/GSTT4
Protein Name:	Glutathione S-transferase theta-1/Glutathione S-transferase theta-4
Other Name:	GSTT1; Glutathione S-transferase theta-1; GST class-theta-1; Glutathione transferase T1-1; GSTT4; Glutathione S-transferase theta-4; GST class-theta-4
Human Swiss Prot No.:	P30711/A8MPT4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GSTT1/4 Polyclonal Antibody detects endogenous levels of GSTT1/4 protein.
Storage and Stability:	-20°C/1 year