

Antibodies Product Datasheet

Product Name:	GRASP55 Antibody
Catalog Number:	C15940
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GORASP2. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	GORASP2
Protein Name:	Golgi reassembly-stacking protein 2
Other Name:	GORASP2; GOLPH6; Golgi reassembly-stacking protein 2; GRS2; Golgi phosphoprotein 6; GOLPH6; Golgi reassembly-stacking protein of 55 kDa; GRASP55; p59
Human Swiss Prot No.:	Q9H8Y8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRASP55 Polyclonal Antibody detects endogenous levels of GRASP55 protein.
Storage and Stability:	-20°C/1 year