

Antibodies Product Datasheet

Product Name:	AGS3 Antibody
Catalog Number:	C15959
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human AGS3.
Concentration:	1 mg/ml
Gene Name:	GPSM1
Protein Name:	G-protein-signaling modulator 1
Other Name:	GPSM1; AGS3; G-protein-signaling modulator 1; Activator of G-protein signaling 3
Human Swiss Prot No.:	Q86YR5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AGS3 Polyclonal Antibody detects endogenous levels of AGS3 protein.
Storage and Stability:	-20°C/1 year