

Antibodies Product Datasheet

Product Name:	Grap Antibody
Catalog Number:	C15967
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRAP. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	GRAP
Protein Name:	GRB2-related adapter protein
Other Name:	GRAP; GRB2-related adapter protein
Human Swiss Prot No.:	Q13588
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Grap Polyclonal Antibody detects endogenous levels of Grap protein.
Storage and Stability:	-20°C/1 year