

Antibodies Product Datasheet

Product Name:	CRIF1 Antibody
Catalog Number:	C15969
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GADD45GIP1. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	GADD45GIP1
Protein Name:	Growth arrest and DNA damage-inducible proteins-interacting protein 1
Other Name:	GADD45GIP1; PLINP1; PRG6; Growth arrest and DNA damage-inducible proteins-interacting protein 1; CKII beta-associating protein; CR6-interacting factor 1; CRIF1; Papillomavirus L2-interacting nuclear protein 1; PLINP; PLINP-1; p53-responsive
Human Swiss Prot No.:	Q8TAE8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRIF1 Polyclonal Antibody detects endogenous levels of CRIF1 protein.
Storage and Stability:	-20°C/1 year