

Antibodies Product Datasheet

Product Name:	GRB2 Antibody
Catalog Number:	C15972
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRB2. AA range:141-190
Concentration:	1 mg/ml
Gene Name:	GRB2
Protein Name:	Growth factor receptor-bound protein 2
Other Name:	GRB2; ASH; Growth factor receptor-bound protein 2; Adapter protein GRB2; Protein Ash; SH2/SH3 adapter GRB2
Human Swiss Prot No.:	P62993
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRB2 Polyclonal Antibody detects endogenous levels of GRB2 protein.
Storage and Stability:	-20°C/1 year