

Antibodies Product Datasheet

Product Name:	GHITM Antibody
Catalog Number:	C15977
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GHITM. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	GHITM
Protein Name:	Growth hormone-inducible transmembrane protein
Other Name:	GHITM; DERP2; MICS1; TMBIM5; My021; Growth hormone-inducible transmembrane protein; Dermal papilla-derived protein 2; Mitochondrial morphology and cristae structure 1; MICS1; Transmembrane BAX inhibitor motif-containing protein 5
Human Swiss Prot No.:	Q9H3K2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GHITM Polyclonal Antibody detects endogenous levels of GHITM protein.
Storage and Stability:	-20°C/1 year