

Antibodies Product Datasheet

Product Name:	GIDRP88 Antibody
Catalog Number:	C15980
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GIDRP88. AA range:721-770
Concentration:	1 mg/ml
Gene Name:	R3HCC1L
Protein Name:	R3H and coiled-coil domain-containing protein 1-like
Other Name:	R3HCC1L; C10orf28; GIDRP88; R3H and coiled-coil domain-containing protein 1-like; Growth inhibition and differentiation-related protein 88; Putative mitochondrial space protein 32.1
Human Swiss Prot No.:	Q7Z5L2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GIDRP88 Polyclonal Antibody detects endogenous levels of GIDRP88 protein.
Storage and Stability:	-20°C/1 year