

Antibodies Product Datasheet

Product Name:	Gas8 Antibody
Catalog Number:	C15994
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Gas8.
Concentration:	1 mg/ml
Gene Name:	GAS8
Protein Name:	Growth arrest-specific protein 8
Other Name:	GAS8; GAS11; Growth arrest-specific protein 8; GAS-8; Growth arrest-specific protein 11; GAS-11
Human Swiss Prot No.:	O95995
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gas8 Polyclonal Antibody detects endogenous levels of Gas8 protein.
Storage and Stability:	-20°C/1 year