

Antibodies Product Datasheet

Product Name:	GFRP Antibody
Catalog Number:	C15997
Reactivity:	Human, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GCHFR. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	GCHFR
Protein Name:	GTP cyclohydrolase 1 feedback regulatory protein
Other Name:	GCHFR; GFRP; GTP cyclohydrolase 1 feedback regulatory protein; GFRP; GTP cyclohydrolase I feedback regulatory protein; p35
Human Swiss Prot No.:	P30047
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GFRP Polyclonal Antibody detects endogenous levels of GFRP protein.
Storage and Stability:	-20°C/1 year