

Antibodies Product Datasheet

Product Name:	GIMAP4 Antibody
Catalog Number:	C16004
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GIMAP4. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	GIMAP4
Protein Name:	GTPase IMAP family member 4
Other Name:	GIMAP4; IAN1; IMAP4; MSTP062; GTPase IMAP family member 4; Immunity-associated nucleotide 1 protein; IAN-1; hIAN1; Immunity-associated protein 4
Human Swiss Prot No.:	Q9NUV9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GIMAP4 Polyclonal Antibody detects endogenous levels of GIMAP4 protein.
Storage and Stability:	-20°C/1 year